

# Sepsis Early Detection Systems (SEDS) User Manual

The Sepsis Early Detection System (SEDS) employs advanced AI predictive models for the diagnosis of sepsis patients, aiming to assist clinicians in diagnosing patients. After successfully registering an account and logging into the system, clinicians can view various basic information of patients, and diagnose patients based on AI model prediction results and medical experience.

Upon logging into SEDS, each clinician will be randomly assigned to a numbered patient. For example, consider patient number 10078. The initial interface clinicians encounter, as shown in Figure 1, provides the related information for the assigned patient. Clinicians are required to complete the rapid bedside diagnosis and treatment plan within 4 hours. Once this is done, the system will display the patient's subsequent examination data, as depicted in Figure 2. This process helps ensure timely and informed decision-making in the treatment of sepsis patients.

Note : The patient's examination information listed below is for reference, intended to illustrate the operation process and system functions. The actual data is subject to the officially deployed system data.

1. 注册

2. 脓毒症

3. AI 模型 -> 脓毒症

4. 预测脓毒症概率0h:1.0,3h:1.0

10078

5. 基础信息

6. 体温: 36.61°C 7. 心率: 88.0bpm

8. 血压: 120/64mmHg

9. 24小时内体液(入量/出量/0ml)

11. 呼吸频率 19insp/min 12. 收缩压: 115mmHg

10. 意识: 意识正常

13. 历史数据

14. 主诉:

15. 中文: 困惑, 呼吸困难 转入重症监护病房的原因: 低血压  
英文: Chief Complaint: Confusion, difficulty breathing Reason for MICU transfer: Hypotension

16. 人口学信息及病史

17. 性别: 男 18. 年龄: 54

19. 体重: 85.4kg 20. 婚姻: 丧偶

21. 种族: 白人 22. 既往病史和现病史

23. 当前时间: 2149/12/11 16:46~2149/12/11 17:46

24. 历史检查数据

25. 初步诊断: 请选择

26. 请选择

27. 4小时内初步治疗方案:

28. 药品: 请选择药品

29. 剂量: 20 盒

30. 药物稀释溶液:

31. 无

32. 500 ml

33. 确定

34. 返回患者目录

35. 下一个患者

Figure 1. Rapid Bedside Diagnosis Interference

1. 下一步检查 [^](#)

2. 降钙素原

3. 血常规

4. 动脉血气分析

5. 止凝血

6. 影像检查

7. 病原检查

8. 培养/涂片

SOFA [^](#)

9. 呼吸系统

Pao2/FiO2:

10. 机械通气:

11. 凝血系统

12. 血小板:

13. 参考范围: 125~350|10<sup>9</sup>/L

14. 肝脏

15. 胆红素: 51.3umol/L

16. 参考范围: 1.71~21|umol/L

17. 心血管系统

MAP: 87.0mmHg

18. 参考范围: 70~105|mmHg

19. 中枢神经系统

gcs: 15.0

20. 参考范围: 3~15

21. 肾脏

22. 肌酐: 97.24umol/L

23. 参考范围: 57~97|umol/L

24. 尿量 (24小时内):

25. 26.

\*最终诊断:

27. 最终4小时内治疗方案:

28. 药品: 29.

30. 用量: 20 31. 盒

32. 药物稀释溶液:

33. 无  500 ml

34. 真实临床上处理相似患者需要花费的时间:

35. 15分钟  36. 5小时

37.

38.  39.

1. Next inspection
2. Procalcitonin
3. Routine blood test
4. Arterial blood gas analysis
5. Anticoagulant
6. Imaging examination
7. Pathogenic examination
8. Culture/Smear
9. Respiratory system Pao2/FiO2
10. Mechanical ventilation
11. Coagulation system
12. Platelet
13. Reference scope : 125~350|10<sup>9</sup>/L
14. Liver
15. Bilirubin : 51.3umol/L
16. Reference scope : 1.71~21|umol/L
17. Cardiovascular system
- MAP : 87.0mmHg
18. Reference scope : 70~105|mmHg
19. Central nervous system
- gcs : 15.0
20. Reference scope : 3~15
21. Kidney
22. Creatinine : 97.24umol/L
23. Reference scope : 57~97|umol/L
24. Urine volume
25. Early diagnosis
26. Please choose
27. Early treatment plan within 4 hours
28. Drugs
29. Please select drugs
30. Dose
31. Box
32. Diluted solutions of drugs
33. None
34. The time required to manage similar patients in real clinical practice
35. Minutes (15 Minutes)
36. Hours (5 Hours)
37. Determine
38. Back to Patients Catalog
39. Next patient

Figure 2. Early Diagnosis Interference

The detailed introduction of the Sepsis Early Detection System (SEDS) is as follows.

# Diagnostic Patient Interface

## Patient ID

After logging into the system account, the clinician randomly assigns a patient with the number 10078. (The patient's examination information listed below is for reference, intended to illustrate the operation process and system functions. The actual data is subject to the officially deployed system data.)



Figure 3. Patient ID

## Model Diagnosis

Diagnosis result of patient model numbered 10078: The patient was diagnosed with sepsis using an AI model. The model predicts a probability of 1 for sepsis to occur currently, and a probability of 1 for sepsis to occur within the next 3 hours.




Figure 4. Model Diagnosis

## Basic Information

The basic information of the patient mainly includes: temperature, blood pressure, heart rate, 24-fluids (input/output), QSOFA score, QSOFA examination items (respiratory rate, consciousness, systolic blood pressure), and historical data of the above examinations. In the historical data, the historical examination results of each examination item are sorted in descending order according to the examination time.

基础信息

体温： 36.61℃      心率： 88.0bpm

血压：                120/64mmHg

24小时内体液(入量/出量):  
量):

QSOFA: 0              意识： 意识正常

呼吸频 19insp/min    收缩压： 115mmHg  
率:

历史数据

Figure 5. Basic Information

历史基础数据

体温 ∨	
2111-10-13 22:36:00	36.94℃
2111-10-13 20:00:00	36.89℃
血压 ∨	
2111-10-13 22:36:00	128/96mmHg
2111-10-13 22:01:00	157/95mmHg
心率 ∨	
2111-10-13 22:36:00	71bpm
2111-10-13 22:00:00	75bpm
输入 ∨	
2111-10-13 21:15:00	氯化钠0.45% NaCl 0.45%/1000.0ml
2111-10-13 20:30:00	葡萄糖50% Dextrose 50%/25.0ml
输出 ∨	
2111-10-13 21:00:00	保险套导尿管 Condom Cath/350.0ml
2111-10-13 18:22:00	保险套导尿管 Condom Cath/85.0ml
QSOFA ∨	
2111-10-13 22:36:00	1
意识 ∨	
2111-10-13 22:36:00	意识改变
2111-10-13 20:00:00	意识正常
呼吸频率 ∨	
2111-10-13 22:36:00	15insp/min
2111-10-13 22:00:00	18insp/min
收缩压 ∨	
2111-10-13 22:36:00	128mmHg
2111-10-13 22:01:00	157mmHg

Figure 6. Historical Data

Demographic Information and Medical History

The demographic information of the patient mainly includes gender, age, weight, marriage, and race. Additionally, click on 'Past Medical History and Present Medical History' to view the patient's medical history.



Figure 7. Demographic Information



Figure 8. Past Medical History and Present Medical History

## Current Time

The current time indicates that the above basic information of the patient is valid during this time. The current time for the patient's basic information is from 16:46 on December 11, 2149 to 17:46 on December 11, 2149.

当前时间:  
2149/12/11 16:46~2149/12/11 17:46

Figure 9. Current Time

## Historical Inspection Data

The historical examination data includes examination items before the current time, mainly including routine blood routine, arterial blood gas analysis, anticoagulation, imaging examination, pathogenic examination, culture/smear, and historical medication. The inspection items in the historical inspection data are based on descending time series, displaying the corresponding inspection item data at historical time points.



Figure 10. Historical Inspection Data

## Rapid Bedside Diagnosis

The clinician makes a rapid bedside diagnosis based on demographic information, medical history, basic data of the current time, and historical data before the current time. The rapid bedside diagnosis mainly confirms whether the patient has sepsis or the suspected degree of sepsis. There are two ways for clinicians to complete rapid bedside diagnosis:

- ① Diagnosis in the selection list: The current system's diagnosis list includes severe sepsis, general sepsis, highly suspected sepsis, low-grade suspected sepsis, and no sepsis. If the clinician's current diagnostic needs are met, the diagnosis name in the list can be directly selected.
- ② Autonomous diagnosis addition: If the diagnosis in the list does not meet the current diagnosis requirements, the clinician can enter the diagnosis name in the blank box below the "Rapid bedside diagnosis" button, click the "Determine" button in the figure, and save the input current diagnosis.
- Note: After completing the treatment plan, you can click the "Determine" button below the "Rapid Bedside Treatment Plan within 4 Hours" to save both the diagnosis and treatment plan.

当前时间:  
2111/10/13 22:36~2111/10/13 23:36

历史检查数据 ▾

\*初步诊断: 请选择

严重脓毒症  
一般脓毒症  
高度疑似脓毒症  
低度疑似脓毒症  
无脓毒症

\*4小时内初步

药品: 请选择

用量: 20

药物稀释溶液

无 ▾ 500 ml +

Figure 11. Rapid Bedside Diagnosis - List Selection

\*初步诊断: 请选择

Figure 12. Rapid Bedside Diagnosis – Self Add

### Rapid Bedside Diagnosis Treatment Plan

Based on the rapid bedside diagnosis results, a rapid bedside diagnosis treatment plan is formulated, which requires filling in the medication, dosage, and drug dilution solution. There are two ways for clinicians to complete the rapid bedside diagnosis treatment plan:

- ① Select List Drug Name: Select the drug name to be used in the drug list, match it by entering a single word in the space box at the top of the list, which helps to efficiently filter the matching drug names in the list. Then select the drug dosage and unit, confirm whether to dilute the solution and dosage, and click the "plus sign" on the right to complete the addition of the drug name to the list.



If the input drug name cannot be matched in the list, use the following method to add the drug name independently.

② Add drug name independently: Enter a custom drug name in the space box at the top of the list without matching it with a drug, then select the drug dosage and unit, confirm whether to dilute the solution and dosage, and finally click the "plus sign" on the right to complete the custom drug name addition.

\*4小时内初步治疗方案:

药品:

请选择药品

无  
抗感染药物-单用  
抗病毒药物  
抗真菌药物  
抗细菌药物 (β-内酰胺类)  
抗细菌药物 (大环内酯类)  
抗细菌药物 (氨基糖甙类)  
抗细菌药物 (喹诺酮类)

Figure 13. Drug List

\*4小时内初步治疗方案:

药品:

用量:

药物稀释溶液:

无

▼

500

ml

+

确定

Figure 14. Autonomous Input Therapy Plan

After the rapid bedside diagnosis and rapid bedside treatment plan within 4 hours, there are two situations in the system operation, which are respectively explained:

- ① If there is a lack of rapid bedside diagnosis or rapid bedside treatment plan within 4 hours, simply click the "Determine" button, and the system will prompt that the current patient's rapid bedside diagnosis or treatment is not completed, and the next steps of the patient's examination cannot be displayed.
- ② After ensuring that the rapid bedside diagnosis and rapid bedside treatment plan within 4 hours are filled in correctly, click the "Determine" button, and the system will record the rapid bedside diagnosis time in the background. At the same time, the system will automatically display the patient's next examination data, and the clinician can continue to diagnose the patient based on

the displayed next examination data.

Special note: Only after completing the rapid bedside diagnosis and the rapid bedside treatment plan within 4 hours and clicking " Determine ", will the system display the patient's next examination data. If the rapid bedside diagnosis or the rapid bedside treatment plan within 4 hours has not been completed, the system will hide the patient's next examination data.



Figure 15. Rapid Bedside Diagnosis Not Filled In



Figure 16. Rapid Bedside Treatment Plan Not Filled In

## Next Step Inspection

After completing the bedside diagnosis and treatment plan within 4 hours, the system will display the patient's next examination data. In the next step of the examination, it mainly includes: procalcitonin, blood routine, arterial blood gas analysis, anticoagulation, imaging examination, pathogen examination, and culture/smear. The clinician further diagnoses the patient's condition by selecting the necessary examination items. It should be noted that the time required for the next examination may affect the timeliness of the patient's treatment.



Figure 17. Next Step Inspection

## SOFA

Clinicians can view the SOFA examination items of the patient during the current time, which mainly include six major parts: respiratory system, coagulation system, liver, cardiovascular system, central nervous system, and kidneys. The respiratory system data needs to be displayed after arterial blood gas analysis, and the coagulation system data needs to be displayed after blood routine examination.

SOFA ^	
呼吸系统	
Pao2/FiO2:	
机械通气:	
凝血系统	
血小板:	
参考范围:	125-350(10 <sup>9</sup> /L)
肝脏	
胆红素:	51.3umol/L
参考范围:	1.71~21umol/L
心血管系统	
MAP:	87.0mmHg
参考范围:	70~105mmHg
中枢神经系统	
gcs:	15.0
参考范围:	3~15
肾脏	
肌酐:	97.24umol/L
参考范围:	57~97umol/L
尿量 (24小时内):	

Figure 18. SOFA

## Early diagnosis

The clinician makes a early diagnosis of the patient based on basic information and all data from the next examination. There are two ways for the clinician to complete the early diagnosis:

- ① Diagnosis in the selection list: The current system's diagnosis list includes severe sepsis, general sepsis, highly suspected sepsis, low-grade suspected sepsis, and no sepsis. If the clinician's current diagnostic needs are met, the diagnosis name in the list can be directly selected.
- ② Self add diagnosis: If the diagnosis in the list does not meet the current diagnostic needs, the clinician can enter the diagnosis name in the blank box below the "Diagnosis" button.

Note: After completing the early 4-hour treatment plan, you can click the " Determine " button to save both the diagnosis and treatment plan.

\*最终诊断: 请选择

请选择

严重脓毒症

一般脓毒症

高度疑似脓毒症

低度疑似脓毒症

无脓毒症

\*最终4小时

药品: 请选择

用量: 20

药物稀释溶液: 无 500 ml +

Figure 19. Early Diagnosis - List Selection

\*最终诊断: 请选择

Figure 20. Early Diagnosis - Self Add

## Early Treatment Plan

After determining the early diagnosis, it is necessary to develop a early treatment plan within 4 hours. There are two ways for clinicians to complete the early treatment plan:

- ① Select List Drug Name: Select the drug name to be used in the drug list, match it by entering a single word in the space box at the top of the list, which helps to efficiently filter the matching drug names in the list. Then select the drug dosage and unit, confirm whether to dilute the solution and dosage, and click the "plus sign" on the right to complete the addition of the drug name to the list. If the input drug name cannot be matched in the list, use the following method to add the drug

name independently.

② Add drug name independently: Enter a custom drug name in the space box at the top of the list without matching it with a drug, then select the drug dosage and unit, confirm whether to dilute the solution and dosage, click the "plus sign" on the right to complete the custom drug name addition.

\*最终4小时内治疗方案:

药品:

用量:

药物稀释溶液:

Figure 21. Early Treatment Plan Within 4 hours

## Time Assessment and Determination

After completing the early diagnosis and treatment plan, clinicians must choose the time required to treat similar patients in real clinical practice. By pulling the progress bar to select the length of time, after completing all the above steps, click "confirm" to complete the diagnosis of the patient! Special note: When clicking "confirm", there may be two situations:

- ① If there is a lack of early diagnosis or treatment plan within 4 hours, simply click the "Determine" button, and the system will prompt that the current patient's early diagnosis or treatment is incomplete, and the patient's early diagnosis and treatment plan data cannot be saved.
- ② Ensure that the bedside diagnosis and early treatment plan are filled in correctly within 4 hours, and complete the time evaluation. Click the "Determine" button, and the system will record the early diagnosis time in the background, successfully saving the patient's early diagnosis, treatment plan, and time evaluation data.

真实临床上处理相似患者需要花费的时间:

15分钟  5小时

Figure 22. Time Assessment and Determination

\*最终诊断: 请选择

\*最终4小时内治疗方案:

药品: 请选择药品

尚未进行诊断!

确定 关闭网页

真实临床上处理相似患者需要花费的时间:

15分钟 5小时

确定

返回患者目录 下一个患者

Figure 23. Early Diagnosis Not Filled In

\*最终诊断: 严重脓毒症

严重脓毒症

\*最终4小时内治疗方案:

药品: 请选择药品

尚未选择治疗方案!

确定 关闭网页

真实临床上处理相似患者需要花费的时间:

15分钟 5小时

确定

返回患者目录 下一个患者

Figure 24. Early Treatment Plan Not Filled In

## **Next Patient**

There are two possible outcomes of clicking on the 'Next Patient' operation, which explain:

- ① Successfully transitioned to the next patient. This situation indicates that there are still undiagnosed patients allocated by the system in the current clinician's account, and the system selects undiagnosed patients for display.
- ② Cannot jump to the next patient. This situation indicates that all patients allocated by the system in the current clinician's account have been diagnosed.